


Integrated computational pipeline for the complete sequencing of glycopeptide using HCD/CID-MS/MS and MS3

Mingqi Liu, Wenfeng Zeng, Yang Zhang, Jianqiang Wu, Simin He, Pengyuan Yang

Institutes of Biomedical Sciences, Fudan University, Shanghai
Institute of Computing Technology, Chinese Academy of Sciences, Beijing

CNCP-2014, 11.12-13

Biological function of site-specific glycosylation: Example of IgG



Long term goal: Interpretation of intact glycopeptide in complex sample

- Study of site-specific glycosylation is developing rapidly
- Identification of intact glycopeptide using LC-MS is the method of choice at the moment
- “Small” obstacle: low abundance, low ionization efficiency
- “Big” obstacle: incomplete information from single MS spectrum

PSM : peptide - spectrum match

“GPSM” : glycopeptide - (2~3) spectra match

Long term goal:
Interpretation of intact glycopeptide in complex sample

- ***Possible strategies***
- Our pipeline
- Challenges




Complete sequencing of glycopeptide: CID + ETD?

- Typical fragmentation methods:


CID: collision based, in trap

HCD: collision based, in beam-type module

ETD: electric based




Glycopeptide spectra : ETD



- ETD is widely used in O-glycopeptide profiling. Large-scale N-glycopeptide profiling using ETD is less common.

- It seems that ETD has the “m/z” barrier around 1,000 ~ 1,100. No such barrier is observed in CID.

Glycopeptide spectra : ETD






m/z :
1415.6

m/z :
1132.7

m/z :
944.1

Glycopeptide spectra : CID




- De-glycosylation is required. Application is limited to N-glycosylation.

- Strictly speaking, information of peptide backbone is incomplete

Long term goal:
Interpretation of intact glycopeptide in complex sample

- Possible strategies
- ***Our pipeline***
- Challenges




Integrated computational pipeline for the complete sequencing of glycopeptide using HCD/CID-MS/MS and MS3





- Sample:
Tryptic glycopeptides from mixture of seven standard glycoproteins

- Instrument:
Thermo LTQ Orbitrap Fusion

HCD-MS/MS @ 40ev : Interpretation of the mass of peptide backbone




HCD-MS/MS @ 40ev : Interpretation of the mass of peptide backbone



- Interpretation of HCD-MS/MS is protein database free



HCD-MS/MS @ 40ev : Interpretation of the mass of peptide backbone




- Target-decoy analysis on HCD-MS/MS. The cutoff we used is 0.4.

CID-MS/MS: Interpretation of the glycan


MS2 #26667 RT: 40.43 AV: 1 NL: 8.42E5
T: FTMS + p NSI d Full ms2 1032.45@cid35.00 [279.00-2000.00]




CID-MS/MS: Interpretation of the glycan



- Glycopeptide CID-MS/MS:
 - 1) Extensive glycosidic bond cleavage
 - 2) Intensity of individual fragments is variable



Interpretation of peptide sequence from MS3





Interpretation of peptide sequence from MS3

- Appropriate target and MS parameter for MS3
- We have tested

CID-MS/MS-DDA-MS3
CID-MS/MS-DIA-MS3
HCD-MS/MS-DDA-MS3
HCD-MS/MS-DIA-MS3

CID or HCD?



Interpretation of peptide sequence from MS3

- HCD-MS/MS based MS3 can only be performed on Orbitrap Fusion at the moment*

Before Orbitrap Fusion:




IT: ion trap
(CID, MS_n)

OT: Orbitrap
(Detector)

Q: Quadrupole
(HCD)

Orbitrap Fusion:



Interpretation of peptide sequence from MS3


- Conventional DDA MS3 from HCD-MS/MS could collect a lot of spectra with peptide sequence information
- Peptide database searching software can be used for peptide sequence in MS3 (pFind, Mascot, SEQUEST, Maxquant, etc.)
- We have develop an integrated software package for the complete sequencing of glycopeptide from MS/MS and MS3

Interpretation of peptide sequence from MS3: Spectral quality


Identified peptide backbone sequence	100 semi	500 semi	2500 semi	All semi
CANLVPVPITJATLDR	Y	Y		
EEQFJSTFR	Y	Y	Y	
EEQFJSTYR	Y	Y		
EEQYJSTYR	Y	Y	Y	
EJGTISR	Y			
EJGTVSR				
FSIJTTNVMGTS�TVR	Y			
FYFTPJK	Y	Y	Y	
FYFTPJKTEDITFLR	Y	Y	Y	Y
GLTFQQJASSM	Y	Y	Y	Y
GNEANYYSJATTDEHGLVQF	Y	Y		
JNSDISSTR	Y	Y		
JSSYLNVR	Y			
LAGKPTHVJVS	Y	Y	Y	Y
LAGKPTHVJVSVVM	Y	Y	Y	Y
LAGKPTHVJVSVVMAEVDGTC	Y			
LAGKPTHVJVSVVMAEVDGTCY	Y	Y		
LVPVPITJATLDQITGK	Y	Y	Y	Y
LVPVPITJATLDR	Y	Y	Y	Y
MJVSVCGLYTYGK	Y			
MVSHHJLTTG	Y	Y	Y	
MVSHHJLTTGA	Y	Y	Y	Y
MVSHHJLTTGAT	Y	Y	Y	Y
MVSHHJLTTGATLINEQWLLTTAK	Y	Y		
NEEYJK	Y	Y	Y	Y
NLFLJH	Y	Y	Y	
NLFLJHSEJATAK	Y	Y	Y	Y

- Spectral quality of MS3 was tested by changing the size of protein database*
- The quality MS3 spectra is sufficient for protein database searching


Complete sequencing of glycopeptide




- For each glycopeptide identification, a pair of HCD/CID-MS/MS is needed
- A group of glycopeptide identification shares one or more MS3 spectra for peptide backbone sequencing




Exemplary spectra: CID-MS/MS



Exemplary spectra: HCD-MS/MS



Exemplary spectra: MS3



Long term goal:
Interpretation of intact glycopeptide in complex sample

- Possible strategies
- Our pipeline
- ***Challenges***

Application of our pipeline to complex sample: Potential challenge

- MS/MS spectral processing

Protein-database free algorithm should work fine

- MS/MS spectral collection

We intentionally split MS/MS and MS3 spectral collection into two MS runs, therefore the balance between quality and quantity should work well

- MS3 spectral processing


Conventional protein database searching of tryptic peptide

- MS3 spectral collection

Is the sensitivity of MS3 enough?



Application of our pipeline to complex sample: Challenge in MS3

- For MS/MS, the latest instrument shows almost zero difference in terms of sensitivity between IT and OT.
- However, for MS3, it is completely different.



*AGC target: maximum
Accumulation time: 500ms*

Application of our pipeline to complex sample:
Top 500 MS3 TIC - IT vs OT



Conclusion

- We have developed Integrated computational pipeline for the complete sequencing of glycopeptide using HCD/CID-MS/MS and MS3
- We will try to interpret glycopeptide in complex sample after some fine tuning in MS3 spectral collection

Acknowledgement

- IBS, Fudan University
Pr. Pengyuan Yang
Yang Zhang
Pan Fang
Guoquan Yan
Weiqian Cao
- ICT, CAS
Pr. Simin He
Wenfeng Zeng
Jianqiang Wu
Hao Chi
Ruixiang Sun
- Thermo Fisher Scientific, Shanghai
Aiyang Nie